# GEM Premier Blood Gas Analyzer and iQM System

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
# Quality Control Guidelines

- CHM.34400 Daily QC - Blood Gas Instruments
  - A minimum of **one level of quality control for pH, pCO<sub>2</sub> and pO<sub>2</sub>** is analyzed at **least every eight hours** of operation when patient specimens are tested, or more frequently if specified in the manufacturer's instructions or laboratory procedure, and when changes occur that may impact patient results.
- CHM.34500 Daily QC - Blood Gas Instruments
  - **The control materials for pH, pCO<sub>2</sub> and pO<sub>2</sub> represent both high and low values on each day of patient testing.**
- CHM.34600 QC - Blood Gas Instruments
  - **At least one level of quality control material for pH, pCO<sub>2</sub> and pO<sub>2</sub> is included each time patient specimens are tested**, except for automated instruments that internally calibrate at least once every 30 minutes of use.
- ISO 15189 section 5.6.1
  - **Design internal quality control systems that verify attainment of the intended quality goals**
- CLSI document EP23 Laboratory QC Based on Risk Management
  - **Establish test method quality requirements based on risk management**

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
([www.cms.hhs.gov/clia/cowppmp.asp](http://www.cms.hhs.gov/clia/cowppmp.asp))

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## Measured Analytes

Analyte	Displayed Ranges
pH	6.80 to 7.80
pCO <sub>2</sub>	5 to 115 mmHg
pO <sub>2</sub>	0 to 760 mmHg
Na <sup>+</sup>	100 to 200 mmol/L
K <sup>+</sup>	0.1 to 20.0 mmol/L
Ca <sup>++</sup>	0.10 to 5.00 mmol/L
Glu	20 to 500 mg/dL
Lac	0.3 to 15.0 mmol/L
Hct	15% to 65%

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	<b>GEM 3500</b>	<b>GEM 4000</b>
<b>Dimensions</b>	<b>44.5x33.0x30.0 cm</b>	<b>45.5x30.5x38.1 cm</b>
<b>Weight</b>	<b>14.2 kg / 1.9 kg</b>	<b>20 kg / 3.6 kg</b>
<b>Analyte Menu</b>	<b>BG, Hct</b>	<b>BG, Hct, Hb, tBili</b>
	<b>BG, Lytes, Hct</b>	<b>BG, Lytes(+ Cl<sup>-</sup>), Hct, Hb, tBili</b>
	<b>BG, Lytes, Glu, Lac, Hct</b>	<b>BG, Lytes(+ Cl<sup>-</sup>), Glu, Lac, Hct, Hb, tBili</b>
<b>Sample Volume</b>	<b>BG/Hct: 135 µL</b>	<b>All item: 150 µL</b>
	<b>BG/Hct/Lytes: 135 µL</b>	<b>CO-Ox: 100 µL</b>
	<b>All item (capillary): 145 µL</b>	<b>All item w/o CO-Ox (Capillary): 65 µL</b>
	<b>All item: 150 µL</b>	
<b>Time To Results</b>	<b>85 sec</b>	<b>70 sec (w/o CO-Ox)</b>
		<b>95 sec (w/ CO-Ox)</b>

1. BG = pH, pCO<sub>2</sub>, pO<sub>2</sub>
2. Lytes = Na<sup>+</sup>, K<sup>+</sup> and Ca<sup>++</sup>
3. Hb = tHb, O<sub>2</sub>Hb, HHb, COHb, MetHb, sO<sub>2</sub>
4. tBili = PAKs available with or without tBili.

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cap

## AQ

# Critical Care/Aqueous Blood Gas With Chemistry

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MENU	PAK SIZE					Onboard Use-life
	75	150	300	450	600†	
Blood Gas, Hct, tHb, O <sub>2</sub> Hb, HHb, Methb, sO <sub>2</sub> , Total Bili* <b>BG</b>	✓	✓	✓	✓	✓	30 days
Blood Gas, Electrolytes, Hct, tHb, O <sub>2</sub> Hb, HHb, Methb, sO <sub>2</sub> , Total Bili* <b>BG+Lyte</b>	✓	✓	✓	✓	✓	30 days
Blood Gas, Electrolytes, Glu, Lac, Hct, tHb, O <sub>2</sub> Hb, HHb, Methb, sO <sub>2</sub> , Total Bili* <b>BG+Lyte+Glu/Lac</b>	✓	✓	✓	✓	✓	30 days
Menus with: Creat, BUN, measured tCO <sub>2</sub> are in development.						


† Onboard use-life is 21 days.

\* Cartridges available with or without Total Bili.

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



- 感應卡與試劑的出廠檢驗
- 卡匣的追溯性 → NIST (National Institute of Standards and Technology)

## Certificate of Traceability and Uncertainty

**Product:** GEM® Premier 4000 Cartridge

**Approach:** per ISO 17511 with uncertainty expressed with K = 2 (95% confidence)

Analytes	Methods and Materials
General approach	Analytes are assigned against working calibrators using validated standing measurement procedures. Working calibrators are assigned using validated selected measurement procedures which are based on the following reference methods:
pH & Ca <sup>++</sup>	per CLSI no. C46-A for pH and CLSI no. C39-A for iCa, i.e., direct potentiometry versus secondary standards prepared from: • phosphate salts for pH, NIST #1861 & 186II, and • calcium carbonate for Ca <sup>++</sup> , NIST #915.
pCO <sub>2</sub> & pO <sub>2</sub>	per CLSI nos. C21-A & C46-A, i.e., tonometry at 37°C using gas mixtures assured to ±.01% versus NIST
Na <sup>+</sup> & K <sup>+</sup>	per NIST pub. no. 260-60 (Aug.'78) and 260-63 (May'79), i.e., flame photometry versus secondary standards prepared from NaCl & KCl, NIST #919 & NIST #918.
Cl <sup>-</sup>	per NIST pub. no. 260-67 (nov.'79) and CLSI no. RS10-P, i.e., coulometric-amperometric titration with silver ion versus secondary standards prepared from NIST #919.
Glucose & Lactate	by automated spectrophotometry: • glucose: hexokinase method per CDC no. 77-8330 and CLSI no. RS1-A using secondary standard prepared from NIST #917 • lactate: no reference method established for lactate. Lactate oxidase is used with secondary standard prepared from USP cat. no. 1614308
tBili	by automated, modified Jendrassik-Grof using standards prepared from NIST #916
Hct	by centrifugation using whole blood per CLSI H7-A3 for establishing correlation. Maintained from lot to lot by controlling conductivity by controlling sodium concentration.
tHb	per CLSI H15-A3 (2000), Hemoglobincyanide by Colorimetry using CyanMetHb standards per CLSI specifications.

### Uncertainty of Measurement at Normal Concentrations

analyte	pH	pCO <sub>2</sub>	pO <sub>2</sub>	Na <sup>+</sup>	K <sup>+</sup>	Ca <sup>++</sup>	Cl <sup>-</sup>	Glucose	Lactate	tHb	tBili	Hct
Units	(eq/kg)	mmHg	mmHg	mmol/L	mmol/L	mmol/L	mmol/L	mg/dL	mmol/L	g/dL	mg/dL	%
U	.005	0.9	1.6	0.8	.06	.011	0.7	1.9	.02	0.2	+0.5	0.5

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

## iQM Features and Technical Description

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All results from 8-hour period require review

VS.

iQM assures quality continuously

iQM reduces error detection time from hours to minutes.<sup>1,2</sup>

	pH	pO <sub>2</sub>	pCO <sub>2</sub>	Na <sup>+</sup>	K <sup>+</sup>	Ca <sup>++</sup>	Glu	Lac	Hct
<b>iQM</b>	<b>3 mins</b>	<b>3 mins</b>	<b>3 mins</b>	<b>17 mins</b>	<b>3 mins</b>	<b>3 mins</b>	<b>11 mins</b>	<b>6 mins</b>	<b>3 mins</b>
Traditional QC (Manual or Auto)	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs	≥ 8 hrs

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HOME "WESTGARD RULES" ESSAYS QC APPLICATIONS LESSONS CLIA & QUALITY DOWNLOADS STORE RES

JAMES WESTGARD  
FOUNDER

Blog

About Us

Reference Materials  
& Resources

CALCULATORS

QC Tools

QC Calculators

HOME ▶ ESSAYS ▶ TRENDS ▶ A MOMENTOUS HAPPENING - A NEW WAY TO DO QC!

## A MOMENTOUS HAPPENING - A NEW WAY TO DO QC!

Written by James O. Westgard

A new QC technology has been cleared that replaces the use of traditional external quality controls. This happened because a manufacturer submitted a claim for this new QC technology and provided the documentation to defend that claim. Read all about it!

***Replaces the use of traditional external QC.***

- The CLIA standard for QC
- A momentous happening!
- What is iQM?
- How does iQM work?
- How does a manufacturer validate QC performance?
- What performance is expected for iQM?
- How does iQM performance compare with traditional QC?
- References

AdChoices ▶

### Espresso ELN

Electronic Lab  
Notebook for  
R&D. Download  
your \*free\* copy  
now!

[www.chembytes.com](http://www.chembytes.com)

(<http://www.westgard.com/essay47.htm>)

20

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*FDA-cleared Intended Use Statement\**

“iQM is an **active** quality process control program designed to provide **continuous monitoring** of the analytical process with **real-time, automatic error detection, automatic correction** of the system and **automatic documentation** of all corrective actions, **replacing the use of traditional external quality controls.**”\*\*

\* Approved in 2002 for GEM Premier 3000 and again in 2006 for GEM Premier 4000.

\*\* Facilities should follow local, state and federal regulatory guidelines to ensure that a total quality management system is followed.

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## Validation of iQM Active Process Control Technology

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



PC Solution	Measurement Assessment	Frequency
A	<ul style="list-style-type: none"> <li>• Used for sensitivity check on all sensors</li> </ul>	Every 4 hours
B	<ul style="list-style-type: none"> <li>• Used to assess sensor sensitivity and performance</li> <li>• Used as corrective actions in high frequency after interference</li> <li>• All sensor outputs are checked every 30 seconds</li> </ul>	After every sample or every 30 minutes
C	<ul style="list-style-type: none"> <li>• Low level of <math>pO_2</math>, pH and <math>pCO_2</math> sensors</li> <li>• Used for conditioning the interference rejection membrane for glucose/lactate sensors</li> </ul>	24 hours
D	<ul style="list-style-type: none"> <li>• Used to assess sensor sensitivity</li> </ul>	12 hours

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
my/PP 03 11/18/2013 10:49 IQM® On 445 Tests 3 Days

**Processing 01:24** Disconnected

Start New Sample Enter Information View Results

JUSTINE, THYME Age: 8 years Mixed Venous, Not Validated

**Micro clot detected for:**  
**Na<sup>+</sup>.**  
**iQM correction in progress.**

OK

**Clot pattern detected.**  
**iQM correction in progress.**

OK

09/05/2002 15:13:24

DATABASE CONFIGURATION CARTRIDGE DIAGNOSTICS SHUTDOWN


pH OK Busy Arterial  
pCO2 OK iQM: On

**Checking for presence of interference and micro clots.**  
**Please wait for results.**

Glucose OK QC  
Lac OK CVP  
Hct OK  
CO-ox OK

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Date/Cart.Lot #	Event	Corrective Action	Result
11/18/2013 11:16:11 0123456789A	Interference Detected After Sample # . Operator:	Operator Notified. Sensor Output Adjusted	Cleared
11/18/2013 11:14:37 0123456789A	Solution A Error for Na <sup>+</sup>	Sensor Output Adjusted	Corrected
11/18/2013 11:14:36 0123456789A	Micro Clot Detected	Sensors Rinsed and Checked	Cleared

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


my/PP 03 11/18/2013 10:49 IQM<sup>®</sup> On 445 Tests 3 Days

Menu Processing 01:24 Disconnected

Start New Sample Enter Information View Results

JUSTINE, THYME Age: 8 years Mixed Venous, Not Validated


**Interference detected for:  
Na<sup>+</sup>, Ca<sup>++</sup>, K<sup>+</sup>  
on last sample, likely due to  
Benzalkonium.**

**! Sensor interference detected for Na+ on last sample.**

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Date/Cart.Lot #	Event	Corrective Action	Result
11/18/2013 11:34:15 0123456789A	Temporary sensor error for Lac during post-sample iQM process.	Operator notified.	Operator acknowledged.
11/18/2013 11:25:26 0123456789A	Interference Detected After Sample # . Operator:	Operator Notified. Sensor Output Adjusted	Cleared
11/18/2013 11:18:31 0123456789A	Solution A Error for Cl	Sensor Output Adjusted	Corrected

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






November 2013      GEM Premier 4000/my/PP 03/888003218

Date/Cart.Lot #	Event	Corrective Action	Result
11/18/2013 11:34:15 0123456789A	Temporary sensor error for Lac during post-sample iQM process.	Operator notified.	Operator acknowledged.
11/18/2013 11:25:26 0123456789A	Interference Detected After Sample # . Operator:	Operator Notified. Sensor Output Adjusted	Cleared
11/18/2013 11:18:31 0123456789A	Solution A Error for Cl <sup>-</sup>	Sensor Output Adjusted	Corrected

page 1 of 7

Print      Close

iQM Corrective Action Report      12/22/2008 10:59:35

Other Months      Print      Exit

Analyzer: GEM 3000, S/N: 12602, Name:      Month: Jul 2002

07/21/2002 23:30:03 Cartridge Lot No.: 213966  
Cartridge Removed.  
Samples Remaining = 412.  
No. of Solution B Adjustments = 6


07/20/2002 07:00:39 Cartridge Lot No.: 213966  
Interference Detected After Sample 38.  
Operator HORMOZ  
Operator Notified. Sensor Output Adjusted.  
Cleared

07/20/2002 06:44:05 Cartridge Lot No.: 213966  
CVP Error for Hct. Operator HORMOZ  
Repeated CVP  
Corrected

07/18/2002 18:29:40 Cartridge Lot No.: 213966  
Cartridge Removed.  
Samples Remaining = 445.  
No. of Solution B Adjustments = 0

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
DATABASE CONFIGURATION CARTRIDGE DIAGNOSTICS SHUTDOWN

1. Select Analytes 2. Type 3. Go

**i** You must run all required levels of CVP material.

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





©2014 Instrumentation Laboratory. All rights reserved. Instrumentation Laboratory is a company of Werfen Group.

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**Accept button**

**Menu CVP Due**

1. Analytes: pH, CT, SIB, pO<sub>2</sub>, pCO<sub>2</sub>, Ca<sup>++</sup>, O<sub>2</sub>Hb, pO<sub>2</sub>, Hct, C/Hb, Na<sup>+</sup>, K<sup>+</sup>, Hct, Hct, K<sup>+</sup>, Lac, TBil

2. Type: CVP /  STAT

3. Go: **GO!**

**Menu CVP Due**


3. Go: **GO!**

Select CVP lot to run or scan the ampoule barcode.

Accepted

- GEM CVP 1 Lot No: 1861 Exp. Date: Refer to Ampoule Label
- GEM CVP 2 Lot No: 2460 Exp. Date: Refer to Ampoule Label
- GEM CVP 3 Lot No: 3854 Exp. Date: Refer to Ampoule Label
- GEM CVP 4 Lot No: 4854 Exp. Date: Refer to Ampoule Label
- GEM CVP 5 Lot No: 5850 Exp. Date: Refer to Ampoule Label

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
The dashboard displays a grid of analyzer cards. Each card shows the analyzer model, location, status, and test/days counts. A search bar is at the top, and a sidebar on the left contains navigation icons for Analyzers, Dashboard, Samples, Operators, Reports, Manage, and Config.

Analyzer Model	Location	Status	Tests	Days
GEM4000 IQM	ICU / G4K-60	Ready	448	13
GEM5000 IQM <sub>12</sub>	ICU3 / G5K-v59	Ready	449	27
NovaNet	Nova / Novanet	Ready	-	-
GEM5000 IQM <sub>12</sub>	East / GEM1	Ready	449	27
GEM5000 IQM <sub>12</sub>	East / GEM2	Ready	448	13
GEM5000 IQM <sub>12</sub>	East / GEM3	Ready	449	27
GEM4000 IQM	Med / G4K	Ready	448	13
GEM5000 IQM <sub>12</sub>	Med G5k	Ready	449	27
GEM4000 IQM	ICU 2 / GEM1	Ready	448	13
GEM5000 IQM <sub>12</sub>	ICU 2 / GEM2	Ready	449	27
GEM5000 IQM <sub>12</sub>	ICU 2 / GEM3	Ready	449	27
GEM4000 IQM	OR1 / GEM 4000	Ready	448	13
GEM5000 IQM <sub>12</sub>	OR2 / GEM 5000	Ready	449	27
GEM4000 IQM	G4KNEW / G4K-60	CVP Due	450	19
GEM5000 IQM <sub>12</sub>	G5K / G5K-62old	-	450	20
GEM5000 IQM <sub>12</sub>	ICU / G5K	Locked	450	20

- Cartridge Status
- Ongoing operations
- Warnings and Errors

**GEM5000 IQM<sub>12</sub>**  
**OR2 / GEM 5000**  
**Ready**  
**449 TESTS** **27 DAYS**

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